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1. Royal Botanic Gardens Kew. 2020. State of the world's plants and fungi. Available at: <https://www.kew.org/science/state-of-the-worlds-plants-and-fungi> [Accessed: 07 December 2021]
2. Finkers, R. Know your onion — The impact of long reads on large genomes. Presentation. Available at: <https://nanoporetech.com/resource-centre/know-your-onion-impact-long-reads-large-genomes> [Accessed: 07 December 2021]
3. United Nations. 2017. World population projected to reach 9.8 billion in 2050, and 11.2 billion in 2100. Available at: <https://www.un.org/development/desa/en/news/population/world-population-prospects-2017.html> [Accessed: 07 December 2021]
4. NCBI National Center for Biotechnology Information. Genomes information by organism. Available at: <https://www.ncbi.nlm.nih.gov/genome/browse#!/overview/> [Accessed: 07 December 2021]
5. Fleischmann, A. et al. Evolution of genome size and chromosome number in the carnivorous plant genus *Genlisea* (*Lentibulariaceae*), with a new estimate of the minimum genome size in angiosperms. *Annals of Botany*. 114 (8): 1651–1663 (2014).
6. Pellicer, J., Fay, M.F., and Leitch, I. J. The largest eukaryotic genome of them all? *Botanical Journal of the Linnean Society*. 164 (1) (2010).
7. Oxford Nanopore Technologies. Ultra-Long DNA Sequencing Kit. Available at: <https://store.nanoporetech.com/ultra-long-dna-sequencing-kit.html> [Accessed: 07 December 2021]
8. Belser, C. et al. Telomere-to-telomere gapless chromosomes of banana using nanopore sequencing. *Commun Biol*. 4(1):1047 (2021).
9. Michael, T. The complex architecture of plant T-DNA transgene insertions. Presentation. Available at: <https://nanoporetech.com/resource-centre/complex-architecture-plant-t-dna-transgene-insertions> [Accessed: 07 December 2021]
10. Long, E.M. et al. Genome-wide imputation using the Practical Haplotype Graph in the heterozygous crop cassava. *bioRxiv* 443913 (2021).
11. Bhat, J.A. et al. Features and applications of haplotypes in crop breeding. *Commun Biol*. 4(1):1266 (2021).
12. Feschotte, C., Jiang, N. and Wessler, S.R. Plant transposable elements: where genetics meets genomics. *Nat Rev Genet* 3(5):329–41 (2002).
13. Schnable, P.S. et al. The B73 maize genome: complexity, diversity, and dynamics. *Science*. 326:1112–1115 (2009).
14. Lee, S.-I and Kim, N.-S. Transposable elements and genome size variations in plants. *Genomics Inform*. 12(3): 87–97 (2014).
15. Debladis, E., Llauro, C., Carpentier, M., Mirouze, M. and Panaud, O. Detection of active transposable elements in *Arabidopsis thaliana* using Oxford Nanopore Sequencing technology. *BMC Genomics*. 18:537 (2017).
16. Jiao, W.B and Schneeberger, K. The impact of third generation genomic technologies on plant genome assembly. *Current Opinion in Plant Biology* 36: 64–70 (2017).
17. Rousseau-Gueutin, M. et al. Long-read assembly of the *Brassica napus* reference genome Darmor-bzh. *Gigascience*. 9(12): giaa137 (2020).
18. Usadel, B. (2017) Complex tomato genomes: Easy with nanopores. Presentation. Available at: <https://nanoporetech.com/index.php/talk/complex-tomato-genomes-easy-nanopores> [Accessed: 07 December 2021]
19. Michael, T.P. et al. High contiguity *Arabidopsis thaliana* genome assembly with a single nanopore flow cell. *Nat. Commun*. 9(1):541 (2018).
20. International Food Policy Research Institute. International germplasm exchanges are getting more complicated, and here's why that matters to global food security. 2017. Available at: <http://www.ifpri.org/blog/international-germplasm-exchanges-are-getting-more-complicated-and-heres-why-matters-global> [Accessed: 07 December 2021]
21. Norris, A.L. et al. Nanopore sequencing detects structural variants in cancer. *Cancer Biol. Ther*. 17(3): 246–253 (2016).
22. Michael, T.P. (2017). Personal communication with Oxford Nanopore Technologies on 29 August 2017.
23. Wittenberg, A. PromethION sequencing of complex plant genomes. Presentation. Available at: <https://nanoporetech.com/resource-centre/PromethION-sequencing-complex-plant-genomes> [Accessed: 07 December 2021]
24. Bolger, A. LOGAN: Lossless graph-based analysis of NGS data sets. Presentation. Available at: <https://nanoporetech.com/resource-centre/logan-lossless-graph-based-analysis-ngs-datasets> [Accessed: 07 December 2021]

25. Steijger, T. et al. Assessment of transcript reconstruction methods for RNA-seq. *Nat. Methods* 10, 1177-1184 (2013).
26. Timp, W. & Jain, M. Direct RNA cDNA sequencing of the human transcriptome. Presentation. Available at: <https://nanoporetech.com/resource-centre/direct-rna-cdna-sequencing-human-transcriptome> [Accessed: 07 December 2021]
27. Ravindran, S. Barbara McClintock and the discovery of jumping genes. *Proc Natl Acad Sci U S A.* 109(50):20198-9 (2012).
28. He, H.-J., Chen, T. and Zhu, J.-K. Regulation and function of DNA methylation in plants and animals. *Cell Res.* 21(3): 442-465 (2011).
29. Gao, Y. et al. Quantitative profiling of N6-methyladenosine at single-base resolution in stem-differentiating xylem of *Populus trichocarpa* using nanopore direct RNA sequencing. *Genome Biol.* 22:22 (2021).
30. Ni, P. et al. Genome-wide detection of cytosine methylations in plant from nanopore data using deep learning. *Nat Commun.* 12(1):5976 (2021).
31. Naish, M. et al. The genetic and epigenetic landscape of the *Arabidopsis* centromeres. *Science.* 374(6569) (2021).
32. Henderson, I. The genetic and epigenetic landscape of the *Arabidopsis* centromeres. Presentation. Available at: <https://nanoporetech.com/resource-centre/video/ncm21/the-genetic-and-epigenetic-landscape-of-the-arabidopsis-centromeres> (2021) [Accessed: 8 December 2021]
33. Pucker, B. Effective characterization of T-DNA insertion lines through nanopore sequencing. Presentation. Available at: <https://nanoporetech.com/resource-centre/video/lc21/effective-characterization-of-t-dna-insertion-lines-through-nanopore-sequencing> (2021) [Accessed: 8 December 2021]
34. Pucker, B., Kleinbölting, N., and Weisshaar, B. Large scale genomic rearrangements in selected *Arabidopsis thaliana* T-DNA lines are caused by T-DNA insertion mutagenesis. *BMC Genomics* 22(599) (2021).
35. NCBI. Col-0_GKat-w. Genome assembly. Available at: https://www.ncbi.nlm.nih.gov/assembly/GCA_905067165.1 [Accessed: 8 December 2021]
36. Choi, J. Y. Unlocking the plant 3D genome architecture with Pore-C sequencing. Presentation. Available at: <https://nanoporetech.com/resource-centre/video/ncm21/unlocking-the-plant-3d-genome-architecture-with-pore-c-sequencing> [Accessed: 7 December 2021]
37. Ulahannan, N. et al. Nanopore sequencing of DNA concatemers reveals higher-order features of chromatin structure. *bioRxiv* 833590 (2019).
38. Oxford Nanopore Technologies. Pore-C: multi-contact, chromosome conformation capture for both genome-wide and targeted analyses. Poster. Available at: <https://nanoporetech.com/resource-centre/porec>. [Accessed: 17 December 2021]
39. Dong et al. Genome-wide Hi-C analysis reveals extensive hierarchical chromatin interactions in rice. *Plant J.* 94(6):1141-1156 (2018).
40. Choi, J. Y. et al. Ancestral polymorphisms shape the adaptive radiation of *Metrosideros* across the Hawaiian Islands. *Proc. Natl. Acad. Sci. U.S.A.* 118(37): e2023801118 (2021).
41. Wittenberg, A. Accuracy improvements in crop genome assembly using the Q20+ chemistry. Available at: <https://nanoporetech.com/resource-centre/video/ncm21/accuracy-improvements-in-crop-genome-assembly-using-the-q20-plus-chemistry> [Accessed: 14 December 2021]
42. KeyGene News. Available at: <https://www.keygene.com/news-events/press-release-keygene-revolutionizes-crop-innovation-using-nanopore-sequencing-technology/> [Accessed: 14 December 2021]
43. Chawla, H.S. Long reads reveal small-scale genome structural variations in allotetraploid canola. Presentation at London Calling 2019.
44. Chawla, H.S. et al. Long-read sequencing reveals widespread intragenic structural variants in a recent allopolyploid crop plant. *Plant Biotechnol J.* 19(2):240-250 (2021).
45. GitHub. NGMLR. Available at: <https://github.com/philres/ngmlr> [Accessed: 07 December 2021]
46. GitHub. Sniffles. Available at: <https://github.com/fritzsedlazeck/Sniffles> [Accessed: 07 December 2021]

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Oxford Nanopore Technologies、风轮图标、EPI2ME、Flongle、GridION、Metrichor、MinION、MinKNOW、PromethION、SmidgION、Ubik和VolTRAX是Oxford Nanopore Technologies plc在不同国家的注册商标。包含的所有其他品牌和名称均为其各自所有者的财产。© 2021 Oxford Nanopore Technologies plc。版权所有。